

For Research Use Only

CLEC3B Recombinant antibody, PBS Only (Capture)

Catalog Number: 85689-1-PBS



Basic Information

Catalog Number:

85689-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003278.2

GeneID (NCBI):

7123

UNIPROT ID:

P05452

Full Name:

C-type lectin domain family 3,
member B

Calculated MW:

23 kDa

Purification Method:

Protein A purification

CloneNo.:

243017D7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02030-1, Tetranectin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85689-1-PBS. Detection antibody: 85689-2-PBS. Standard: Eg3063. Range: 0.313-40 ng/mL