

For Research Use Only

HOMER1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85692-2-PBS



Basic Information

Catalog Number:

85692-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3109

GenBank Accession Number:

BC015502

GeneID (NCBI):

9456

UNIPROT ID:

Q86YM7

Full Name:

homer homolog 1 (Drosophila)

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

243007B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02028-1, HOMER1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85692-3-PBS. Detection antibody: 85692-2-PBS. Standard: Ag3109. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP02028-2, HOMER1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85692-1-PBS. Detection antibody: 85692-2-PBS. Standard: Ag3109. Range: 0.781-100 ng/mL