

For Research Use Only

HOMER1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85692-5-PBS



Basic Information

Catalog Number:

85692-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3109

GenBank Accession Number:

BC015502

GeneID (NCBI):

9456

UNIPROT ID:

Q86YM7

Full Name:

homer homolog 1 (Drosophila)

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

243007B4

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02028-3, Human HOMER1 Recombinant Matched Antibody Pair - PBS only. 85692-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag3109. 85692-5-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean HOMER1 concentration was determined to be 6.1 ng/mL in LNCaP cell extract based on a 2.6 mg/mL extract load.