

For Research Use Only

DLGAP3 Recombinant antibody

Catalog Number: 85695-1-RR



Basic Information

Catalog Number:

85695-1-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001080418

GeneID (NCBI):

58512

UNIPROT ID:

O95886

Full Name:

discs, large (Drosophila) homolog-associated protein 3

Calculated MW:

106 kDa

Observed MW:

106-120 kDa

Purification Method:

Protein A purification

CloneNo.:

243052B5

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Background Information

DLGAP3, also named as DAP3, belongs to the SAPAP family. It may play a role in the molecular organization of synapses and neuronal cell signaling. DLGAP3 could be an adapter protein linking ion channel to the subsynaptic cytoskeleton. It may induce enrichment of PSD-95/SAP90 at the plasma membrane.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85695-1-RR (DLGAP3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.