

For Research Use Only

RLBP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85697-4-PBS



Basic Information

Catalog Number:

85697-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7602

GenBank Accession Number:

BC004199

GeneID (NCBI):

6017

UNIPROT ID:

P12271

Full Name:

retinaldehyde binding protein 1

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

242989A1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02104-1, RLBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85697-4-PBS. Detection antibody: 85697-3-PBS. Standard: Ag7602. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP02104-2, RLBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85697-4-PBS. Detection antibody: 85697-2-PBS. Standard: Ag7602. Range: 0.625-80 ng/mL