

For Research Use Only

DIS3 Recombinant antibody, PBS Only

Catalog Number: 85705-2-PBS



Basic Information

Catalog Number:

85705-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31725

GenBank Accession Number:

NM_014953

GeneID (NCBI):

22894

UNIPROT ID:

Q9Y2L1

Full Name:

DIS3 mitotic control homolog (S. cerevisiae)

Calculated MW:

109KD

Observed MW:

110 kDa

Purification Method:

Protein A purification

CloneNo.:

242934G8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85705-2-RR (DIS3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85705-2-PBS in a different storage buffer formulation.