

For Research Use Only

DIS3 Recombinant antibody

Catalog Number: 85705-2-RR



Basic Information

Catalog Number:

85705-2-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31725

GenBank Accession Number:

NM_014953

GeneID (NCBI):

22894

UNIPROT ID:

Q9Y2L1

Full Name:

DIS3 mitotic control homolog (S. cerevisiae)

Calculated MW:

109KD

Observed MW:

110 kDa

Purification Method:

Protein A purification

CloneNo.:

242934G8

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HEK-293T cells, Jurkat cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85705-2-RR (DIS3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.