

For Research Use Only

SPATA2 Recombinant antibody

Catalog Number: 85712-1-RR



Basic Information

Catalog Number:

85712-1-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30773

GenBank Accession Number:

BC009481

GeneID (NCBI):

9825

UNIPROT ID:

Q9UM82

Full Name:

spermatogenesis associated 2

Calculated MW:

58 kDa

Observed MW:

58-60 kDa

Purification Method:

Protein A purification

CloneNo.:

242523H4

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HeLa cells, MCF-7 cells, Raji cells, PC-3 cells, A549 cells, mouse testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85712-1-RR (SPATA2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 85712-1-RR against Human SPATA2 were performed. The affinity constant is 0.175 nM.