

For Research Use Only

# FNIP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85714-3-PBS



## Basic Information

Catalog Number: 85714-3-PBS	GenBank Accession Number: BC001956	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 96459	CloneNo.: 242941C8
Source: Rabbit	UNIPROT ID: Q8TF40	
Isotype: IgG	Full Name: folliculin interacting protein 1	
Immunogen Catalog Number: AG28976	Calculated MW: 508 aa, 57 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP02105-1, FNIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85714-3-PBS. Detection antibody: 85714-2-PBS. Standard: Ag28976. Range: 0.625-80 ng/mL