

For Research Use Only

# CYP2C8 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85717-2-PBS



## Basic Information

Catalog Number: 85717-2-PBS	GenBank Accession Number: BC020596	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 1558	CloneNo.: 242948H6
Source: Rabbit	UNIPROT ID: P10632	
Isotype: IgG	Full Name: cytochrome P450, family 2, subfamily C, polypeptide 8	
Immunogen Catalog Number: AG9898	Calculated MW: 490 aa, 56 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP02106-1, CYP2C8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85717-2-PBS. Detection antibody: 85717-3-PBS. Standard: Ag9898. Range: 1.563-200 ng/mL.