

For Research Use Only

CYP2C8 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85717-4-PBS



Basic Information

Catalog Number:

85717-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9898

GenBank Accession Number:

BC020596

GeneID (NCBI):

1558

UNIPROT ID:

P10632

Full Name:

cytochrome P450, family 2, subfamily C, polypeptide 8

Calculated MW:

490 aa, 56 kDa

Purification Method:

Protein A purification

CloneNo.:

242948D5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02106-2, CYP2C8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85717-4-PBS. Detection antibody: 85717-3-PBS. Standard: Ag9898. Range: 1.563-200 ng/mL.