

For Research Use Only

# C1QA Recombinant antibody, PBS Only (Detector)

Catalog Number: 85719-2-PBS



## Basic Information

Catalog Number:

85719-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_015991.4

GeneID (NCBI):

712

UNIPROT ID:

P02745

Full Name:

complement component 1, q  
subcomponent, A chain

Calculated MW:

26kDa

Purification Method:

Protein A purification

CloneNo.:

250093G7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP02098-2, C1QA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85719-3-PBS. Detection antibody: 85719-2-PBS. Standard: Eg2396. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP02098-1, C1QA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85719-1-PBS. Detection antibody: 85719-2-PBS. Standard: Eg2396. Range: 1.563-200 ng/mL.