

For Research Use Only

ZCCHC3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85721-1-PBS



Basic Information

Catalog Number:

85721-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30047

GenBank Accession Number:

NM_033089

GeneID (NCBI):

85364

UNIPROT ID:

Q9NUD5

Full Name:

zinc finger, CCHC domain containing

3

Calculated MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

250088A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02096-1, ZCCHC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85721-3-PBS. Detection antibody: 85721-1-PBS. Standard: Ag30047. Range: 0.469-60 ng/mL.