

For Research Use Only

# ZCCHC3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85721-2-PBS



## Basic Information

Catalog Number:

85721-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30047

GenBank Accession Number:

NM\_033089

GeneID (NCBI):

85364

UNIPROT ID:

Q9NUD5

Full Name:

zinc finger, CCHC domain containing

3

Calculated MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

250088D1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP02096-2, ZCCHC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85721-3-PBS. Detection antibody: 85721-2-PBS. Standard: Ag30047. Range: 0.469-60 ng/mL.