

For Research Use Only

ZCCHC3 Recombinant antibody

Catalog Number: 85721-3-RR



Basic Information

Catalog Number:

85721-3-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30047

GenBank Accession Number:

NM_033089

GeneID (NCBI):

85364

UNIPROT ID:

Q9NUD5

Full Name:

zinc finger, CCHC domain containing

3

Calculated MW:

44 kDa

Observed MW:

43-45 kDa

Purification Method:

Protein A purification

CloneNo.:

250088G1

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, K-562 cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85721-3-RR (ZCCHC3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 85721-3-RR against Human ZCCHC3 were performed. The affinity constant is below 1 pM.