

For Research Use Only

Mouse Bche Recombinant antibody, PBS Only (Capture)

Catalog Number: 85722-3-PBS



Basic Information

Catalog Number:

85722-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2899

GenBank Accession Number:

NM_009738.3

GeneID (NCBI):

12038

UNIPROT ID:

Q03311

Full Name:

butyrylcholinesterase

Calculated MW:

68 kDa

Purification Method:

Protein A purification

CloneNo.:

250086C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02095-1, MOUSE Bcbe Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85722-3-PBS. Detection antibody: 85722-2-PBS. Standard: Eg2899. Range: 0.625-80 ng/mL.