

For Research Use Only

ALDH3A1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85725-4-PBS



Basic Information

Catalog Number:

85725-4-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7977

GenBank Accession Number:

BC004370

GeneID (NCBI):

218

UNIPROT ID:

P30838

Full Name:

aldehyde dehydrogenase 3 family, memberA1

Calculated MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

243134D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP02087-1, ALDH3A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85725-4-PBS. Detection antibody: 85725-3-PBS. Standard: Ag7977. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP02087-2, ALDH3A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85725-4-PBS. Detection antibody: 85725-2-PBS. Standard: Ag7977. Range: 1.563-200 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 85725-4-PBS against Human ALDH3A1 was performed. The affinity constant is below 1 pM.