

For Research Use Only

ZEB1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85729-3-PBS



Basic Information

Catalog Number:

85729-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16204

GenBank Accession Number:

BC112392

GeneID (NCBI):

6935

UNIPROT ID:

P37275

Full Name:

zinc finger E-box binding homeobox 1

Calculated MW:

1125 aa, 124 kDa

Purification Method:

Protein A purification

CloneNo.:

243139B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02086-2, ZEB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85729-4-PBS. Detection antibody: 85729-3-PBS. Standard: Ag16204. Range: 0.625-80 ng/mL