

For Research Use Only

# Mouse Nov Recombinant antibody, PBS Only (Capture)

Catalog Number: 85731-2-PBS



## Basic Information

Catalog Number:

85731-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2856

GenBank Accession Number:

NM\_010930.4

GeneID (NCBI):

18133

UNIPROT ID:

Q64299

Full Name:

nephroblastoma overexpressed gene

Calculated MW:

39 kDa

Purification Method:

Protein A purification

CloneNo.:

250054H6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02088-1, MOUSE Nov/CCN3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85731-2-PBS. Detection antibody: 85731-1-PBS. Standard: Eg2856. Range: 1.563-200 ng/mL.