

For Research Use Only

DNALI1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85734-4-PBS



Basic Information

Catalog Number:

85734-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11652

GenBank Accession Number:

BC046117

GeneID (NCBI):

7802

UNIPROT ID:

O14645

Full Name:

dynein, axonemal, light intermediate chain 1

Calculated MW:

258 aa, 30 kDa

Purification Method:

Protein A purification

CloneNo.:

243157A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02073-1, DNALI1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85734-4-PBS. Detection antibody: 85734-3-PBS. Standard: Ag11652. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP02073-2, DNALI1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85734-4-PBS. Detection antibody: 85734-2-PBS. Standard: Ag11652. Range: 0.781-100 ng/mL.