

For Research Use Only

C1S Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85736-4-PBS



Basic Information

Catalog Number:

85736-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5977

GenBank Accession Number:

BC056903

GeneID (NCBI):

716

UNIPROT ID:

P09871

Full Name:

complement component 1, s
subcomponent

Calculated MW:

77 kDa

Purification Method:

Protein A purification

CloneNo.:

250031F9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02077-1, C1S Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85736-4-PBS. Detection antibody: 85736-3-PBS. Standard: Ag5977. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP02077-2, C1S Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85736-2-PBS. Detection antibody: 85736-4-PBS. Standard: Ag5977. Range: 0.781-100 ng/mL.