

For Research Use Only

Mouse Dcn Recombinant antibody, PBS Only (Capture)

Catalog Number: 85737-1-PBS



Basic Information

Catalog Number:

85737-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2955

GenBank Accession Number:

NM_001190451.2

GeneID (NCBI):

13179

UNIPROT ID:

P28654

Full Name:

decorin

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

250099G6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02099-2, MOUSE Decorin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85737-1-PBS. Detection antibody: 85737-2-PBS. Standard: Eg2955. Range: 1.563-200 ng/mL.