

For Research Use Only

Mouse PF4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85738-2-PBS



Basic Information

Catalog Number:

85738-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3345

GenBank Accession Number:

NM_019932

GeneID (NCBI):

56744

UNIPROT ID:

Q9Z126

Full Name:

platelet factor 4

Calculated MW:

11kd

Purification Method:

Protein A purification

CloneNo.:

242741B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02069-2, MOUSE CXCL4/PF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 98397-2-PBS. Detection antibody: 85738-2-PBS. Standard: Eg3345. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP02069-1, MOUSE CXCL4/PF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85738-1-PBS. Detection antibody: 85738-2-PBS. Standard: Eg3345. Range: 0.781-100 ng/mL