

For Research Use Only

# HNRNPUL2-Specific Recombinant antibody

Catalog Number: 85739-2-RR



## Basic Information

Catalog Number:

85739-2-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001079559

GeneID (NCBI):

221092

UNIPROT ID:

Q1KMD3

Full Name:

heterogeneous nuclear  
ribonucleoprotein U-like 2

Calculated MW:

85 kDa

Observed MW:

110-120 kDa

Purification Method:

Protein A purification

CloneNo.:

250029E1

Recommended Dilutions:

WB 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : MCF-7 cells, HEK-293T cells, Jurkat cells, HeLa cells, SW 1990 cells, NIH/3T3 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85739-2-RR (HNRNPUL2-Specific antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.