

For Research Use Only

NRG4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85742-2-PBS



Basic Information

Catalog Number:

85742-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0997

GenBank Accession Number:

NM_138573.4

GeneID (NCBI):

145957

UNIPROT ID:

Q8WWG1

Full Name:

neuregulin 4

Calculated MW:

13kDa

Purification Method:

Protein A purification

CloneNo.:

242431A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02113-2, NRG4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85742-3-PBS. Detection antibody: 85742-2-PBS. Standard: Eg0997. Range: 1.25-80 ng/mL.



Cytometric bead array standard curve of MP02113-1, NRG4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85742-1-PBS. Detection antibody: 85742-2-PBS. Standard: Eg0997. Range: 0.625-80 ng/mL.