

For Research Use Only

# TRUB1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85746-2-PBS



## Basic Information

Catalog Number: 85746-2-PBS	GenBank Accession Number: BC030601	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 142940	CloneNo.: 243129G2
Source: Rabbit	UNIPROT ID: Q8WWH5	
Isotype: IgG	Full Name: TruB pseudouridine (psi) synthase homolog 1 (E. coli)	
Immunogen Catalog Number: AG3213	Calculated MW: 349 aa, 37 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02110-1, TRUB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85746-3-PBS. Detection antibody: 85746-2-PBS. Standard: Ag3213. Range: 0.781-100 ng/mL.