

For Research Use Only

TRUB1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85746-3-PBS



Basic Information

Catalog Number: 85746-3-PBS	GenBank Accession Number: BC030601	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 142940	CloneNo.: 243129F8
Source: Rabbit	UNIPROT ID: Q8WWH5	
Isotype: IgG	Full Name: TruB pseudouridine (psi) synthase homolog 1 (E. coli)	
Immunogen Catalog Number: AG3213	Calculated MW: 349 aa, 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02110-1, TRUB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85746-3-PBS. Detection antibody: 85746-2-PBS. Standard: Ag3213. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP02110-2, TRUB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85746-3-PBS. Detection antibody: 85746-1-PBS. Standard: Ag3213. Range: 0.781-100 ng/mL