

For Research Use Only

CARD9 Recombinant monoclonal antibody, PBS Only

Catalog Number: 85749-1-PBS



Basic Information

Catalog Number:

85749-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33643

GenBank Accession Number:

BC008877

GeneID (NCBI):

64170

UNIPROT ID:

Q9H257

Full Name:

caspase recruitment domain family, member 9

Calculated MW:

62 kDa

Observed MW:

60-65 kDa

Purification Method:

Protein A purification

CloneNo.:

250045G3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

CARD9, a member of the caspase recruitment domain (CARD) protein family, plays an important regulatory role in cell apoptosis. Defects in CARD9 are the cause of familial candidiasis type 2 (CANDF2), characterized by impaired clearance of fungal infections and resulting in colonization and infections of the mucosa or skin, predominantly with *Candida albicans*. CARD9 is highly expressed in the spleen and regulates activation of NF- κ B by interaction with the positive regulator of apoptosis BCL10 via the CARD domain.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85749-1-RR (CARD9 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85749-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85749-1-RR against Human CARD9 were performed. The affinity constant is 14.3 nM.