

For Research Use Only

# MATK Recombinant antibody, PBS Only (Detector)

Catalog Number: 85755-2-PBS



## Basic Information

Catalog Number:

85755-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0128

GenBank Accession Number:

BC000114

GeneID (NCBI):

4145

UNIPROT ID:

P42679

Full Name:

megakaryocyte-associated tyrosine kinase

Calculated MW:

56 kDa

Purification Method:

Protein A purification

CloneNo.:

250071F2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02089-2, MATK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85755-3-PBS. Detection antibody: 85755-2-PBS. Standard: Ag0128. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP02089-1, MATK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85755-1-PBS. Detection antibody: 85755-2-PBS. Standard: Ag0128. Range: 1.563-200 ng/mL.