

For Research Use Only

TMOD3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85756-3-PBS



Basic Information

Catalog Number:

85756-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3215

GenBank Accession Number:

BC020542

GeneID (NCBI):

29766

UNIPROT ID:

Q9NYL9

Full Name:

tropomodulin 3 (ubiquitous)

Calculated MW:

352 aa, 40 kDa

Purification Method:

Protein A purification

CloneNo.:

250041C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02090-1, TMOD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85756-3-PBS. Detection antibody: 85756-2-PBS. Standard: Ag3215. Range: 3.125-200 ng/mL.



Cytometric bead array standard curve of MP02090-2, TMOD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85756-3-PBS. Detection antibody: 85756-1-PBS. Standard: Ag3215. Range: 3.125-200 ng/mL.