

For Research Use Only

TMOD3 Recombinant antibody, PBS Only

Catalog Number: 85756-5-PBS



Basic Information

Catalog Number: 85756-5-PBS	GenBank Accession Number: BC020542	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 29766	CloneNo.: 250041D6
Source: Rabbit	UNIPROT ID: Q9NYL9	
Isotype: IgG	Full Name: tropomodulin 3 (ubiquitous)	
Immunogen Catalog Number: AG3215	Calculated MW: 352 aa, 40 kDa	
	Observed MW: 40 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Tropomodulin 3 (TMOD3) belongs to the tropomodulin family, which controls cytoskeletal dynamics by regulating the end-capping of actin filaments. TMOD3 bridges the gap between structural integrity and signaling cascades. TMOD3 regulates cytoskeletal architecture and is implicated in cell adhesion, migration, and mechanotransduction.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85756-5-RR (TMOD3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85756-5-PBS in a different storage buffer formulation.