

For Research Use Only

TNFAIP8L2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85759-1-PBS



Basic Information

Catalog Number:

85759-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7964

GenBank Accession Number:

BC063014

GeneID (NCBI):

79626

UNIPROT ID:

Q6P589

Full Name:

tumor necrosis factor, alpha-induced protein 8-like 2

Calculated MW:

184 aa, 21 kDa

Purification Method:

Protein A purification

CloneNo.:

250047E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02092-1, TNFAIP8L2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85759-1-PBS. Detection antibody: 85759-3-PBS. Standard: Ag7964. Range: 0.781-100 ng/mL.