

For Research Use Only

# RIPK2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85760-6-PBS



## Basic Information

Catalog Number:

85760-6-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0909

GenBank Accession Number:

BC004553

GeneID (NCBI):

8767

UNIPROT ID:

O43353

Full Name:

receptor-interacting serine-threonine kinase 2

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

250069A4

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP02093-3, Human RIPK2 Recombinant Matched Antibody Pair - PBS only. 85760-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0909. 85760-5-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL.

The mean RIPK2 concentration was determined to be 76.00 ng/mL in THP-1 cell extract based on a 1.50 mg/mL extract load.