

For Research Use Only

# Rat VEGFR1/Flt-1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85761-1-PBS



## Basic Information

Catalog Number:	85761-1-PBS	GenBank Accession Number:	Purification Method:
Concentration:	1 mg/ml	GeneID (NCBI):	Protein A purification
Source:	Rabbit	54251	CloneNo.:
Isotype:	IgG	UNIPROT ID:	242311D7
Immunogen Catalog Number:	EG1167	P53767	
		Full Name:	
		fms-related tyrosine kinase 1	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
rat

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MPO2054-1, RAT VEGFR-1/Flt-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85761-3-PBS. Detection antibody: 85761-1-PBS. Standard: Eg1167. Range: 1.563-200 ng/mL.