

For Research Use Only

# DDX23 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85763-3-PBS



## Basic Information

Catalog Number:

85763-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0258

GenBank Accession Number:

BC002366

GeneID (NCBI):

9416

UNIPROT ID:

Q9BUQ8

Full Name:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 23

Calculated MW:

96 kDa

Purification Method:

Protein A purification

CloneNo.:

243147D4

## Applications

Tested Applications:

IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP02081-2, DDX23 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85763-3-PBS. Detection antibody: 85763-2-PBS. Standard: Ag0258. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP02081-1, DDX23 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85763-3-PBS. Detection antibody: 85763-1-PBS. Standard: Ag0258. Range: 0.781-100 ng/mL.