

For Research Use Only

DGKA Recombinant antibody, PBS Only (Detector)

Catalog Number: 85764-2-PBS



Basic Information

Catalog Number:

85764-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2075

GenBank Accession Number:

BC023523

GeneID (NCBI):

1606

UNIPROT ID:

P23743

Full Name:

diacylglycerol kinase, alpha 80kDa

Calculated MW:

735 aa, 83 kDa

Purification Method:

Protein A purification

CloneNo.:

243133G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02079-2, DGKA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85764-3-PBS. Detection antibody: 85764-2-PBS. Standard: Ag2075. Range: 1.563-200 ng/mL.