

For Research Use Only

# PRL Recombinant antibody, PBS Only (Detector)

Catalog Number: 85765-2-PBS



## Basic Information

Catalog Number:	85765-2-PBS	GenBank Accession Number:	BC015850	Purification Method:	Protein A purification
Concentration:	1 mg/ml	GeneID (NCBI):	5617	CloneNo.:	242898B12
Source:	Rabbit	UNIPROT ID:	P01236		
Isotype:	IgG	Full Name:	prolactin		
Immunogen Catalog Number:	AG9764	Calculated MW:	227 aa, 26 kDa		

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP02083-1, Prolactin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85765-3-PBS. Detection antibody: 85765-2-PBS. Standard: Ag9764. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP02083-2, Prolactin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85765-1-PBS. Detection antibody: 85765-2-PBS. Standard: Ag9764. Range: 0.781-100 ng/mL