

For Research Use Only

PRL Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 85765-5-PBS



Basic Information

Catalog Number:

85765-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5868

GenBank Accession Number:

NM_000948.5

GeneID (NCBI):

5617

UNIPROT ID:

P01236

Full Name:

prolactin

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

252191D12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Prolactin concentration was determined to be 8.9 ng/mL in human serum.



Cytometric bead array standard curve of MP02083-3, Prolactin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85765-5-PBS. Detection antibody: 85765-4-PBS. Standard: Eg5868. Range: 0.156-20 ng/mL.