

For Research Use Only

# Cleaved Caspase 4/P20 Recombinant monoclonal antibody

Catalog Number: 85770-1-RR



## Basic Information

Catalog Number:

85770-1-RR

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC017839

GeneID (NCBI):

837

UNIPROT ID:

P49662

Full Name:

caspase 4, apoptosis-related cysteine peptidase

Calculated MW:

377 aa, 43 kDa

Observed MW:

20 kDa

Purification Method:

Protein A purification

CloneNo.:

243150D6

Recommended Dilutions:

WB: 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Staurosporine treated Jurkat cells,

## Background Information

Caspase-4 participates in cleavage and activates gasdermin D (GSDMD), leading to cell membrane perforation, cell swelling and rupture, thereby inducing pyroptosis.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

