

For Research Use Only

# TOMM40 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85773-1-PBS



## Basic Information

Catalog Number:

85773-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13065

GenBank Accession Number:

BC017224

GeneID (NCBI):

10452

UNIPROT ID:

O96008

Full Name:

translocase of outer mitochondrial  
membrane 40 homolog (yeast)

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

250010F3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP02065-2, TOMM40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85773-3-PBS. Detection antibody: 85773-1-PBS. Standard: Ag13065. Range: 0.781-100 ng/mL