

For Research Use Only

PDCD6/ALG2 Recombinant antibody, PBS Only

Catalog Number: 85792-3-PBS



Basic Information

Catalog Number:

85792-3-PBS

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2949

GenBank Accession Number:

BC012384

GeneID (NCBI):

10016

UNIPROT ID:

O75340

Full Name:

programmed cell death 6

Calculated MW:

191 aa, 22 kDa

Purification Method:

Protein A purification

CloneNo.:

250087C5

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 85792-3-PBS against Human PDCD6/ALG2 were performed. The affinity constant is 3.63 nM.