

For Research Use Only

Mouse IL-27A p28 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85795-3-PBS



Basic Information

Catalog Number:

85795-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0897

GenBank Accession Number:

GeneID (NCBI):

246779

UNIPROT ID:

O35228(IL27B)+Q8K3I6(IL27A)

Full Name:

interleukin 27

Purification Method:

Protein A purification

CloneNo.:

242655D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02119-1, MOUSE IL-27 p28/IL-30 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85795-3-PBS. Detection antibody: 85795-2-PBS. Standard: Eg0897. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP02119-2, MOUSE IL-27 p28/IL-30 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85795-3-PBS. Detection antibody: 85795-1-PBS. Standard: Eg0897. Range: 0.781-100 ng/mL