

For Research Use Only

# SPR Recombinant antibody, PBS Only (Detector)

Catalog Number: 85796-1-PBS



## Basic Information

Catalog Number:

85796-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10395

GenBank Accession Number:

BC017310

GeneID (NCBI):

6697

UNIPROT ID:

P35270

Full Name:

sepiapterin reductase (7,8-dihydrobiopterin:NADP+ oxidoreductase)

Calculated MW:

261 aa, 28 kDa

Purification Method:

Protein A purification

CloneNo.:

243138D5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02100-1, SPR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85796-3-PBS. Detection antibody: 85796-1-PBS. Standard: Ag10395. Range: 0.781-100 ng/mL.