

For Research Use Only

# Mouse Cd74 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85802-3-PBS



## Basic Information

Catalog Number:

85802-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001042605.1

GeneID (NCBI):

16149

UNIPROT ID:

P04441-1

Full Name:

CD74 antigen (invariant polypeptide  
of major histocompatibility complex,  
class II antigen-associated)

Calculated MW:

32 kDa

Purification Method:

Protein A purification

CloneNo.:

250109D7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP02137-2, MOUSE CD74 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85802-3-PBS. Detection antibody: 85802-2-PBS. Standard: Eg2716. Range: 0.781-50 ng/mL.