

For Research Use Only

LTBP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85820-2-PBS



Basic Information

Catalog Number:

85820-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25392

GenBank Accession Number:

BC130289

GeneID (NCBI):

4052

UNIPROT ID:

Q14766

Full Name:

latent transforming growth factor beta
binding protein 1

Calculated MW:

1721 aa, 187 kDa

Purification Method:

Protein A purification

CloneNo.:

250127A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02129-2, LTBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85820-2-PBS. Detection antibody: 85820-1-PBS. Standard: Ag25392. Range: 0.625-40 ng/mL