

For Research Use Only

Mouse Tnc Recombinant antibody, PBS Only (Capture)

Catalog Number:85821-3-PBS



Basic Information

Catalog Number: 85821-3-PBS	GenBank Accession Number: N/A	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 21923	CloneNo.: 242866B10
Source: Rabbit	UNIPROT ID: Q80YX1-1	
Isotype: IgG	Full Name: tenascin C	
	Calculated MW: 232 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02127-1, MOUSE TNC/Tenascin-C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85821-3-PBS. Detection antibody: 85821-2-PBS. Standard: Eg2849. Range: 0.313-40 ng/mL.