

For Research Use Only

Mouse uPAR/CD87 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85824-3-PBS



Basic Information

Catalog Number:

85824-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3156

GenBank Accession Number:

NM_011113.4

GeneID (NCBI):

18793

UNIPROT ID:

P35456

Full Name:

plasminogen activator, urokinase receptor

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

242976C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02130-1, MOUSE uPAR/CD87 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85824-4-PBS. Detection antibody: 85824-3-PBS. Standard: Eg3156. Range: 0.625-80 ng/mL