

For Research Use Only

Mouse Ctf1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85828-3-PBS



Basic Information

Catalog Number:

85828-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2867

GenBank Accession Number:

NM_007795.1

GeneID (NCBI):

13019

UNIPROT ID:

Q60753

Full Name:

cardiotrophin 1

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

242890G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02133-1, MOUSE Ctf1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85828-3-PBS. Detection antibody: 85828-2-PBS. Standard: Eg2867. Range: 0.781-25 ng/mL.



Cytometric bead array standard curve of MP02133-2, MOUSE Ctf1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85828-3-PBS. Detection antibody: 85828-1-PBS. Standard: Eg2867. Range: 0.781-25 ng/mL.