

For Research Use Only

C6orf15 Recombinant antibody

Catalog Number: 85830-1-RR



Basic Information

Catalog Number:

85830-1-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20276

GenBank Accession Number:

BC153157

GeneID (NCBI):

29113

UNIPROT ID:

Q6UXA7

Full Name:

chromosome 6 open reading frame 15

Calculated MW:

325 aa, 34 kDa

Observed MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

250126F4

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HaCaT cells, Caco-2 cells, PC-3 cells

Background Information

C6orf15, also named as Protein STG, is a 325 amino acid protein, which as a secreted protein is expressed in skin and tonsils. The calculated molecular weight of C6orf15 is 34 kDa, but we detected a 50-55 kDa protein by western blot. C6orf15 may contains some modified residues.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85830-1-RR (C6orf15 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.