

For Research Use Only

C6orf15 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85830-3-PBS



Basic Information

Catalog Number:

85830-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20276

GenBank Accession Number:

BC153157

GeneID (NCBI):

29113

UNIPROT ID:

Q6UXA7

Full Name:

chromosome 6 open reading frame 15

Calculated MW:

325 aa, 34 kDa

Purification Method:

Protein A purification

CloneNo.:

250126D9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02146-1, C6orf15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85830-4-PBS. Detection antibody: 85830-3-PBS. Standard: Ag20276. Range: 0.781-100 ng/mL.